
# A-Series

Under vehicle lift

SV	5	Installationsanvisning / Bruksanvisning
EN	25	Installation instructions / User manual
DE	45	Einbauanleitung / Benutzerhandbuch
FR	65	Manuel d'installation / Manuel d'utilisation
ES	85	Instrucciones de instalación / Modo de empleo
IT	105	Istruzioni per l'installazione / Manuale delle istruzioni

SV

EN

DE

FR

ES

IT

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Thank you for choosing a

# A-Series from BraunAbility!

The manual that follows is an integral, important part of the product, providing you with information on how to perform accurate and safe installation of the product. Please read through the manual carefully before starting your work.

*Indicates a warning for a situation in which damage, body injury or death could result if the warning is ignored.*

If you have any questions about your customer's adaptation, please feel free to contact us. Your feedback is important to us. It helps us in our work of continuously improving our product and its installation. You can find more information about us and our range of car adaptation solutions by visiting [www.braunability.eu](http://www.braunability.eu).

Once again, thank you for placing your confidence in our products!

*This product is designed for a specific purpose and should only be used accordingly.*

## Installation

- 
- 28 Pallet contents on delivery
  - 29 Safety
  - 30 Preparing the cassette for installation
  - 31 Installing the cassette
  - 32 Adjusting the outer end position
  - 33 Install the pump module
  - 34 Routing hose and cables
  - 36 Routing and connecting the electrical supply
  - 38 Adjusting the upper limit
  - 39 Adjusting the inclination
  - 40 Verification by the installer

## Service check list

---

42

## Technical specification

---

45

## User instruction


---

46


Diagrams and exploded views are at the end of the manual.

## Installation


### Pallet contents on delivery

x1 Lift cassette

x1 Pump module with cover

x1 Pump module handle

x1 Hand-held control


x1 Holder


x2 Clips


x5 Screws

x1 Long red cable


x1 Short red cable


x1 Black cable


x1 Extra cable lug

x1 Plastic protector

x1 Fuse holder


x1 Fuse

x2 Screws

## Safety

- DANGER**  This lift has only been tested and approved for transferring people with reduced mobility into a vehicle.
- DANGER**  The lift must not be used as a link-bridge to the pavement or raised groundwork. In order to achieve optimum embarkation and alighting, the lift platform should be flat on the ground.
- DANGER**  Wheelchair brakes should always be applied during the lifting operation. The front roll stop and the inner barrier, see illustration below, should be raised during operation and be clear of any interference throughout. This is especially important to consider when an assistant is on the platform.
- DANGER**  A safety valve is triggered when the maximum weight capacity of 350 kg is exceeded, when this is activated the platform will not lift. Should this happen, reduce the load and try again.
- DANGER**  A battery switch lock prevents unauthorized use. Turn to OUT to remove the switch.


## Preparing the cassette for installation


*Always use appropriate personal protective equipment during installation.*

- 1 Raise the vehicle to facilitate access during the entire installation. Start by mounting the brackets, with the accompanying fasteners, delivered in a separate package. Follow the instructions supplied with the brackets. Note that the brackets should only be tightened lightly so there is some free play for adjusting the positioning of the cassette.
- 2 Cut away the banding around the cardboard packaging. Remove the cardboard hood, the boxes and the wooden transport brackets. Leave the cassette on the pallet.
- 3 Temporarily remove the screws and associated washers from the four mounting slides located along the sides of the cassette.

## Installing the cassette


- 1 Raise the pallet using a forklift. Position the cassette between the brackets, at least 10 mm from the rear axle.
- 2 Adjust the mounting slides so they correspond with the respective brackets. Fasten the cassette to the brackets using the screws and washers removed in an earlier step. Use the rearmost two holes. When attaching a footstep, an accessory that can be ordered separately, use the same mounting slide. The third hole is used only in this instance.
- 3 Lower the pallet and remove it. Retighten the screws to 50 Nm. Then retighten all the bracket fastenings onto the vehicle chassis.

## Adjusting the outer end position


If the cassette lift is installed in combination with a footstep the outer end position of the lift needs to be adjusted. This to ensure that the bridge plate reaches sufficiently far into the vehicle. Adjust the outer end position before fitting the footrest.


## Install the pump module

- 1 Install the pump module in the desired location inside the vehicle. Remember to place the module on the nearside of the vehicle to facilitate easier manual operation in the event of a lift failure. Use a suitable fastening in each of the three mounting points on the floor and in at least one of the mounting points on the side panel. Ensure the installation is sufficiently stable to withstand operation in an emergency when using the manual hand pump.
- 2 Store the handle mentioned above close to the pump module so it is easily accessible when required. Use the accompanying clips.
- 3 Remove the yellow transport plug from the oil reservoir air hose and replace it with the supplied vent cap.

## Routing hose and cables


- 1 Hydraulic hoses and electrical couplings between the power pack and the cassette, can be routed through the same hole. Use existing holes and spaces in the vehicle chassis for this purpose if possible. If this is not possible, drill a hole in a suitable position with a 50 mm hole saw. Place the plastic grommet into the hole to protect the cables against sharp edges.
- 2 Lay the oil hose from the pump module to the cassette in an appropriate manner. Note that bends in the hose must not be less than 30 mm in diameter in order to ensure correct function. Connect the hose to the quick-coupling on the cassette.

- 3 Lay the wiring loom from the cassette to the pump module in an appropriate manner. The connectors are marked A and B. Plug into the corresponding sockets on the pump module.
- 4 When the lift is fully inserted into the cassette, wire C36 gives a + (plus) 12V signal. If a - (minus) 12V signal is desired via the vehicle ignition, see the connection diagram at the end of this manual.  
Note! Cable C36 is only intended for connecting the lamp.
- 5 Connection to the Interlock kit (item number 451139), or another system that prevents unauthorised use. Separate instructions are included with the kit.

*According to the EN 1756-2:2004+A1:2009 standard, the lift must be equipped with a system that prevents unauthorised use in the absence of the operator. Other local requirements may apply.*


## Routing and connecting the electrical supply


- 1 Fit the fuse holder along with fuse in close proximity to the vehicle battery.
- 2 Connect the long red cable with the angled end lug to the battery switch in the pump module. Route the cable towards the battery through the drilled hole in the vehicle floor and connect it to the fuse holder and make connection secure. The length of the long red cable can be shortened if too long. In this case, use the supplied cable lug that fits on the fuse holder.

Connect the short red cable. Mount one end to the fuse and the other end up to the 12V positive terminal on the vehicle battery.

- 3 Connect the black ground cable. Fasten one end to the pump module and the other end to a suitable grounding point. A suitable grounding point is an area of clean metal on the vehicle chassis, free from dirt and paint.
- 4 Connect the hand-held control to the corresponding socket on the pump module. Fasten the cable with the p clip to minimise strain on the connection. Fit the holder for the hand-held control in close proximity to the lift so that it is easily accessible when needed.

See supplementary illustration on next page.


When electric supply is connected a complete test run of the lift's full movement cycle must be performed. The lift will stop automatically once it has reached its set limits.

*Keep a close eye on the movements of the lift during the cycle test and make adjustments according to the following pages if needed. Neglecting to do this could result in physical damage on vehicle and/or lift.*


## Adjusting the upper limit

- 1 The inner platform (bridge plate) should be level with the floor. If necessary, adjust the upper limit as illustrated below.

## Adjusting the inclination

- 1 The inner platform (bridge plate) should be level with the floor. If necessary, adjust the inclination as illustrated below.
- 2 Adjust the cam. Ensure the wheel is straight in relation to the cam. Failure to do so can cause damage the bridge plate.
- 3 Manually deploy the travelling frame, see separate emergency instruction, in order to ensure the platform moves freely in and out of the cassette. If it does not, it may be because the mounting brackets are over tightened and may have deformed the sides of the cassette. The brackets should be adjusted if this is the case.

## Verification by the installer

According the directive EN 1756-2:2004+A1:2009, the installer shall confirm the compatibility between the lift and the vehicle.

Static test.....	<p><b>Deformation</b></p> <ul style="list-style-type: none"> <li>- The unladen platform is positioned mid-way between ground level and vehicle floor level and measurements are taken of the height of the platform and its angular attitude relative to the vehicle floor.</li> <li>- An evenly spread load of 125 % (438 kg) of the maximum load (350 kg) is applied to the platform, and subsequently removed.</li> <li>- By repeating measurements of the height and attitude of the platform, verify that no permanent deformation has occurred in any part of the tail lift or its attachment to the vehicle which could affect the function of the tail lift.</li> </ul> <p><b>Drift</b></p> <ul style="list-style-type: none"> <li>- An evenly spread load of 125 % (438 kg) of the maximum load is applied to the platform, positioned at floor vehicle floor level. Measurements are taken of the height of the platform and its angular attitude relative to the vehicle floor. These measurements are repeated after a 15 min test period.</li> <li>- Verify that the vertical drift of the platform between the two measurements has not exceeded 15 mm.</li> <li>- Verify that the angular drift of the platform between the two measurements has not exceeded 2°.</li> </ul>
Dynamic test.....	<ul style="list-style-type: none"> <li>- With the maximum load applied to the platform, verify that the tail lift is able to operate throughout its full range of normal lifting, lowering and any tilting movements. If the lift is unable to lift with this weight, the lifting capacity must be increased. Proceed as follows (see also the illustration on the next page):           <ol style="list-style-type: none"> <li>1. Loosen the locking nut. Use a 14 mm (or 9/16") socket or swivel joint wrench.</li> <li>2. Turn the adjusting screw in a clockwise direction. Use a 3/16" Allen key.</li> <li>3. Retighten the locking nut. Use the Allen key to hold the adjusting screw in position when tightening the locking nut.</li> </ol> </li> </ul>
Test of operation and safety functions .....	<ul style="list-style-type: none"> <li>- All functions of the tail lift and operations of all safety devices are verified after the static and dynamic tests have been completed. These tests do not apply to pipe break valves nor non-resettable safety devices such as electrical fuses (These items are the subject of a manufacturer's type test).</li> </ul>

Test to verify that the tail lift cannot lift excessive load .....

- A load of 438 kg (125 % of the maximum load) is applied to the platform, positioned at ground level. Actuate the "UP" control and verify that the platform does not lift (tilt is permissible). If the lift is able to lift with this weight, the lifting capacity must be decreased. Proceed as follows (see also the illustration below):

1. Undo the locking nut. Use a 14 mm (or 9/16") socket or swivel joint wrench.
2. Turn the adjusting screw in a counterclockwise direction. Use a 3/16" Allen key.
3. Retighten the locking nut. Use the Allen key to hold the adjusting screw in position when tightening the locking nut.


Lift model \_\_\_\_\_

Serial number \_\_\_\_\_

Perform date \_\_\_\_\_

Name of dealer \_\_\_\_\_

Signatur, installation engineer \_\_\_\_\_

## Service check list

To ensure maximum service life and functionality please follow the check list below.

Lift model _____	Dealer _____
Serial number _____	Resp. mechanic _____
Vehicle reg. number _____	Owner _____
Vehicle make/model _____	Owner address _____

\*Hydraulic oil, BraunAbility art. no. 440612



\*\*OKS 3751, BraunAbility art. no. 400878








\*\*\*Light oil

<b>EVERY DAY</b>	<b>A</b>	Check that the bridge plate ( <b>A</b> ) goes up and down in the manner intended.
	<b>B</b>	Check that the hand rails with protective plates ( <b>B</b> ) are functioning in the manner intended. Protective plates are an option.
	<b>C</b>	Check that the outer roll stop ( <b>C</b> ) goes up and down in the manner intended.
	<b>D</b>	Check that the cassette ( <b>D</b> ) is free of external damage.
	<b>E</b>	Check that the handle ( <b>E</b> ) for the pump module is in place.
	<b>F</b>	Check that the hand control ( <b>F</b> ) and cable are intact and functioning.
		Check that there is no visible leaking of hydraulic oil.
		Check that instructions are available and in a legible condition.
<b>EVERY MONTH</b>	<b>G</b>	Check that the level of the *hydraulic oil is correct. When the lift is fully inside the cassette, the level must lie between the "stow level" mark and up to a maximum of 5 mm below the mark ( <b>G</b> ).
	<b>H</b>	Check that the fuses ( <b>H</b> ) are intact.
		Check that visible cables and tubing are free of external damage.
		Check that everything works in the manner intended.
<b>EVERY SIX MONTHS</b>		Check that there is no leaking of hydraulic oil or wear and tear of moveable components behind the lower cover of the cassette. Remove the lower cover to perform the inspection.
	<b>I</b>	Check that the two hinges of the hand rails are ***lubricated ( <b>I</b> ).
	<b>J</b>	Check that the four hinges, which connect the foldable parts of the platform, are ***lubricated ( <b>J</b> ).
	<b>K</b>	Check that the three bushings, holding the linkage arm ( <b>K</b> ) to the bridge plate, are ***lubricated.
	<b>L</b>	Check that all four support rollers are not deformed or damaged.
	<b>M</b>	Clean and **lubricate the tracks ( <b>M</b> ).
	<b>N</b>	Clean and ***lubricate the tracks ( <b>N</b> ). Reinstall the lower cover.
	<b>O</b>	Check that the two hinges ( <b>O</b> ) of the front hatch are ***lubricated and that the hatch is closing properly.
		Check that all bracket bolts are properly tightened.
		Check that the emergency release functions.
		Check that bolts, pins and joints for wear/damage and that they are properly secured.

EN

Service check list


# Technical specification

EN

Technical specification


## Product data

Cassette length.....	1220 mm
Cassette width.....	1050 mm
Cassette height.....	140 mm
Usable platform length .....	1210-1340 mm
Usable platform width.....	800 mm
Minimum overlap, vehicle floor .....	20 mm
Minimum width, door opening.....	820 mm
Cassette weight.....	130 kg
Noise level during use.....	less than 70 dB
Pump module length.....	245 mm
Pump module width .....	120 mm
Pump module height.....	453 mm

## Product performance


Moving, vertical cycle .....	1075 mm
Moving from upper edge of cassette and up.....	430 mm
Minimum move up.....	130 mm
Moving from upper edge of cassette and down.....	645 mm
Load capacity .....	350 kg

## User instruction


## Battery switch lock

Take out the switch lock to prevent unauthorised use.

## Manual emergency operation

### Unstow and down

- 1 Open the valve for 3 seconds. Close the valve again.
- 2 Fold down the front cover. Move the wing screw to hold the front cover down.
- 3 Unlock the red catch. Beware of the strike plate that springs out.
- 4 Pull out the platform and unfold.
- 5 Open the valve to lower the platform to the ground. Close the valve again.

### Up and stow

- 6 Work the pump to raise the platform until the arms align with the cassette. Perform step 4, 3 and 2 in reverse order to stow the lift back into the cassette.
- 7 If not working, open the valve and work the pump to remove air from the system. Close the valve and try step 6 again.

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© 2024 BraunAbility Europe AB
BraunAbility Europe AB is certified
BraunAbility Europe AB reserves the right to make alterations without previous notice.